

Z-LINK™

Achieve greater efficiencies with Z-LINK™. Z-LINK™ is Zetta's own remote connectivity tool that allows instant access to your CT, MRI and PET/CT equipment. Our proactive approach allows us to quickly resolve many system errors and issues without the need of onsite service. This can be especially beneficial if your site is in a remote location.



Z-LINK™ 24x7 connectivity allows us to:

- Examine and Diagnose Network Configurations
- View and Analyze Image Artifacts
- View Error Logs and Analyze Data
- Run Diagnostics
- Remote Fix
- System Configuration Changes
- Add and Remove MRI Coils from your Coil List

Magnet Monitoring Including:

- Helium Levels
- Magnet Pressure
- Coldhead Efficiency
- Chiller Temperature
- Shield Temperature
- Water flow
- Water Temperature

Z-LINK™ is HIPAA Compliant with strong passwords, SSL & 256-bit Rijndael AES encryption and provides the highest security level on all remote connectivity sessions from login to logout.